



END OF OPAL-RAD PACS SWITCH TO: iQ-SYSTEM PACS

Konica Minolta will discontinue support for OPAL-RAD PACS. For those looking for a reliable and future-proof solution, now is the perfect time to upgrade to the powerful iQ-SYSTEM PACS. Discover the numerous benefits of the PACS solution from iQ IMAGE.



COST EFFECTIVE

- Patient portal/physician portal incl. remote order entry
- Multi-specialty use (ortho, cardio, radiology, vet)
- Multiple integration options (DICOM, HL7, text based)



DICOM AND DATA

- Easy DICOM tag modification
- DICOM modality worklist used for device planning, study reconciliation or scanning documents



WORKFLOW

- 3-level Hanging Protocol (specifically for radiologists)
- Intelligent previous study management (prefetching and/or automatic query/retrieval) about 2 times faster



IMPORT AND EXPORT

- QR code for patient portal access
- Automatic email notification to referring physician or manual email to request second opinion
- PDF to DICOM and DICOM to PDF converter



IMAGING

- Dedicated orthopedic measurements, e.g. 1-click Gonstead analysis, Cobb, spinal curves and labeling
- Dedicated radiology measurements, e.g. perpendicular and interior angle, ratio, ROI window
- Dedicated cardio measurements, e.g. EF by Simpson, mitral regurgitation
- Dedicated veterinary measurements, e.g. TPLO, TTA, VHS and 3D measurements for TTTG
- Additional image filters, e.g. smoothen, sharpen, remove X-ray grid (mainly for vets)



INTEGRATION

- Easy integration with almost any RIS using GDT, HL7 or other tools
- Easy integration of almost any 3rd party software with iQ-VIEW using iQ-LAUNCHER



Need help? Contact us at:
sales.us@iq-image.com